

- ⑥ Disconnect used drainage bag and immediately connect new sterile bag to the catheter or sheath
- ⑦ Avoid touching the ends of any connecting area
- ⑧ Use the new protective cap to cover the end of the used drainage bag
- ⑨ Empty the used drainage bag and dispose of according to local policies and guidelines
- ⑩ Remove gloves and discard as normal clinical waste
- ⑪ Always use an aseptic and non-touch technique
- ⑫ Use on a single patient during single procedure









Prosys[®] Non Latex Gloves

- Open Prosys pack to reveal one pair of sterile gloves
- Take care not to touch the inner contents
- Decontaminate hands before and after procedure by washing and drying thoroughly
- Follow the principles of an aseptic and non-touch technique
- Avoid touching the ends of any connecting areas
- Always adhere to local policies and guidelines



If the package is damaged, the device must not be used. This device is not suitable for re-sterilisation.

SYMBOL EXPLANATION

 Do not re-use	STERILE EO Sterilised using ethylene oxide	REF Catalogue number
 Do not re-sterilise	 Caution	LOT Batch code
 Consult instructions for use	 Contains or presence of Phthalates	 Date of manufacture
 Do not use if package is damaged	This device is free from natural latex	 Use by date



CLINISUPPLIES Ltd, First Floor, 1 Blackmoor Lane, Croxley Park, Hertfordshire, WD18 8GA, UK
 Tel: +44 (0)20 8863 4168 Fax: +44 (0)20 8426 0768 Email: info@clinisupplies.co.uk
www.clinisupplies.co.uk

EC REP Advena Ltd, Tower Business Centre, 2nd floor,
 Tower Street, Swatar, BKR 4013, Malta.

12/20 v3



Instructions for use

- ① Open Prosys sterile pack taking care not to touch the inner contents
- ② The sterile gloves will now be on top of the drainage bag
- ③ Wash and dry hands thoroughly
- ④ Place on and wear the sterile gloves
- ⑤ Position and prepare the new sterile bag. Do not remove the protective cap at this stage